

CeleMag Clean-up Bead



CeleMag Clean-up Bead utilizes unique magnetic bead-based chemistry for the clean-up of nucleic acids, providing a simple, flexible and highly reproducible protocols for purification and size selection of nucleic acids.

Benefits


- Highly optimized with Celeemics NGS Target Enrichment Kits
- Industry-leading purification yield efficiency
- Perfect alternative to bead-based purification workflows
- Adjustable and consistent size selection
- Automated liquid handler compatible

Applications


- Next generation sequencing
- Sanger sequencing
- Nucleic acid purification

CeleMag Clean-up Bead provides equivalent NGS Library preparation recovery efficiency in comparison to A Bead.

Library Amount After Pre-PCR (ng)




1st Conc. After Pre-PCR (ng/uL)



- Libraries were prepared with the Celeemics Library Prep kit from 100ng of gDNA sheared by Covaris (17.3ng/uL,239bp)
- Measured Recovery Amount with Qubit and size distribution was identified by TapeStation after 1st Pre-PCR.

CeleMag Clean-up Bead provides highly comparable performance to A bead in size selection workflows, achieving consistent DNA size distributions and exactly desired library sizes.



[Purification Protocols_1.8X]

- A. Mix sample with the Beads.
- B. Incubate the solution for 5 minutes.
- C. After the solution is clear, remove the supernatant.
- D. Ethanol wash and elution.



[Size Selection Protocols]

- 1st Beads Selection
- A. Mix sample with Beads
 - B. Incubate the solution for 5minutes.
- 2nd Beads Selection
- A. Mix the supernatant from 1st selection with Beads.
 - B. Incubate the solution for 5 minutes.
 - C. After the solution is clear, remove the supernatant.
 - D. Ethanol wash and elution.

Ordering Information

Catalog #	Description	Volume
CMCB05	CeleMag Clean-up Bead, 5ml	5ml
CMCB50	CeleMag Clean-up Bead, 50ml	50ml
CMCB500	CeleMag Clean-up Bead, 500ml	500ml